



Thank you for your purchase of the VH Dissector. Following are several helpful steps that will be of benefit and that will minimize your start-up time.

If questions arise at any time, feel free to e-mail us at vhdissessor@toltech.net or call 720-859-4140. Please describe the problem as best you can in order for us to be more efficient in our response.

HELPFUL STEPS

1. The CD or DVD installs on your computer and does not have to remain in the drive for future use. Installation time for the CD is less than 5 minutes, the DVD takes approximately 25-30 minutes. A Wait Message may not appear but let the installation process continue.
2. After installation, a **VH Dissector** icon may appear on your desktop, depending upon how your computer is configured. If an icon does not appear, follow these instructions:
 - a. Click on the **START** menu, then on **ALL PROGRAMS**.
 - b. Select the **ToLTech** folder.
 - c. Right click on the **VH Dissector** (Pro or Lite) icon and select **Send to Desktop (Create Shortcut)**.
3. Launch the **VH Dissector**. A Message will ask if you wish to take the **Tutorial** (this message only appears the first time that you access the VH Dissector. If you answer no, the Tutorial may be accessed at a later time by selecting the **Lessons** Tab, then selecting **Tutorial**).
4. **Tutorial**. The Tutorial takes approximately one hour to complete. It is advised that every new user of the VH Dissector take the time to go through the Tutorial. Questions that arise during normal usage of the VH Dissector may be sent to vhdissessor@toltech.net, but you are urged to first complete the Tutorial, where many of the questions will be answered.
 - a. **Red Words** in the Tutorial link to VH Dissector states and will cause the dissection image or the cross section image to change (or both),
 - b. **Blue Links** take you to a different section of the Tutorial for purposes of showing a screen shot or a more detailed description. To return, right click and press the **Back** button.
5. The **Quick Reference Guide** and the **Exercise** are electronic versions of those sent with your VH Dissector CD's and/or DVD's. The documents will help you to become more knowledgeable with the VH Dissector's navigational techniques **AFTER you have taken the TUTORIAL**. The documents may be printed or electronically distributed wherever needed.

Once again, thank you for becoming a customer of Touch of Life Technologies. We look forward to assisting you as you use the VH Dissector.

Customer Service
Touch of Life Technologies
vhdissessor@toltech.net
720-859-4140